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Massage mode enhancing the quality of relaxation

SR Sound sleep, Relaxation Mode

Debussy's 'moonlight'+ sound of bugs and birds at night

It features the scenarios inspired by Debussy's 'moonlight' synonymous with lyrical classical music, and sound of bugs and birds making sound quietly at night. This is a mode that promotes recovery from fatigue and sleep with acupressure of reflexology point related to body relaxation to relieve tension.

SE Stretch Exercise Mode

Mozart's 'Turkish March' + sound of calm waves

The scenarios inspired by Mozart's dynamic 'Turkey March', coupled with the sound of soft waves that disappear with bubbles, help you relieve stress and relax your muscles. This mode provides massage to reflexology point of musculoskeletal system, such as neck, shoulder, arm, waist and knee, like stretching whole body.

BA Brain Activation Mode

'Morning scene' in Peer Gynt's Morning Mood music + babbling sound of streams

The scenario inspired by the morning scene in Peer Gynt's Morning Mood music describing the slowly rising dawn light and constant babbling sound of streams help enhance concentration and creativity. This mode provides massages to reflexology point related to brain and organ memory.

SL Slim Line Mode

Chrysler's joy of love' + water dripping sound

The scenarios inspired by Chrysler's 'joy of love' distinguished for bright and happy mood and sound of refreshing water dripping sound is provided. This is a mode massaging the reflexology point which helps manage the cheek health and swelling.

GF Golf Mode

Bach's 'Goldberg Variations' + sound from grassland of the North Island, New Zealand

The scenarios inspired by Bach's 'Goldberg Variations', a classy symphony loved by music enthusiasts, and the recorded sound from grassland of the North Island, New Zealand make you feel as if you were in the field. This is a mode that massages reflexology point to help golfer improve performance and develop flexibility in the shoulders, arms, legs, waist, etc.

MF Mindfulness Mode

Scenery sound of temple

The scenario helps increase focus on the senses and landscape sound makes you feel as you were in temple. This is a mode that helps reduce stress and brain fatigue by enhancing the effects of meditation that brings stability to the mind.

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Warranty Information

- This product passed rigorous quality control standards of our company.
- The warranty period is 1 year. Free repair will be provided for failures occurring within warranty period.
- However, repair will be charged in the event of failure caused due to negligence of consumer or careless storage.

Warranty Conditions

- If the serial number on the product is not clearly legible, proper service may be difficult.
- The following cases are not covered under warranty repair.
 - If the product was repaired by undesignated company
 - Disassembled/repared by customer
 - Failure caused by falling
 - Failure occurred due to careless storage
 - Failure caused due to natural disaster (typhoon, earthquake, flood, etc.)
 - If there is no warranty attached

※ Specific regulations are subject to change without separate notice because warranty service will be continuously upgraded.

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How to Use

■ Power On

Press and hold the Power/Start button to be ready for operation with the voice guidance saying "Power on".

■ Power Off

When the power is on, press and hold the Power/Start button to turn the power off with the voice guidance saying "Power off".

■ Operation

Press the Power/Start button shortly to activate IAN with the voice guidance saying "Start Acupressure".

■ Stop/Start Operation

Press the Power/Start button shortly during acupressure to stop acupressure with the voice guidance saying "End Acupressure". Pressing the Power button shortly again will start acupressure all over again.

■ Adjustment of Size(Head Size)

Press the V(down) button of Size button to adjust the headfit down to fit the head size of user. Press the ^ (release) button to discharge the air pressure supporting the headfit and release the adjustment. After release, headfit can be adjusted to fit the head size by pressing the V(down) button again. Adjust the Head Fit after stopping operation because Head Fit does not work during acupressure.

■ Acupressure Pin Position Adjustment

Each scenario begins with positioning acupressure. If all areas of earlobe, ear upper part, middle ear(inside of the ear), and outer circumference of the ear felt evenly pressed, the position is correct. If it does not fit well, stop operation, adjust the top and bottom with the Size button and move IAN back and forth into correct position.

■ Voice Guidance, White Noise Level Control

Press the Volume button to adjust the volume of voice guidance. Minimizing the volume will mute the sound. Press the Volume button to adjust the white noise level during the acupressure. Minimizing the volume will mute the sound.

■ Remote app provided

Connecting the IAN with a remote app for your smartphone via Bluetooth makes it even more convenient. Please search for 'IAN remote' in the Play Store (App Store) and install it. Please refer to the app description in the Play Store (App Store) for detailed usage.

■ Mode Selection

Pressing the Mode button will make the following mode(function) selected one after another, along with voice guidance.

Mode	Description	Mode	Description	Mode	Description
SR	Sound sleep,Relaxation Mode	SE	Stretch Exercise Mode	GF	Golf Mode
BA	Brain Activation Mode	SL	Slim Line Mode	MF	Mindfulness Mode
MN	Manual mode, voice guidance on operations of IAN				

■ Scenario Selection

Three scenarios are provided for each mode. Pressing the Scenario button will make the following scenarios selected one after another, along with voice guidance. Adapt to the stepwise scenario and then gradually increase the strength.

Scenario	Description
1	Running time about 5 minutes. Ambient sound after voice guidance
2	Running time about 10 minutes. Ambient sound
3	Running time about 10 minutes. Acupressure strength of main reflexology point is high. Ambient sound

Precautions

Be familiarized with precautions of users and use IAN in a safe and correct way.

When using this product, do not wear glasses, jewelry such as earrings, hairpins, or hearing aids. Please put your hair behind your ears before use.

Individuals with sensitive ears may feel pain in their ears. If fever or pain persists, stop using this product.

Depending on the individual, healing reaction(a temporary worsening of symptoms that occurs when the body is going through the process of healing itself through the process of finding the balance of the body) may occur. Please refer to the healing reaction.

Do not use the product with wet hands or face.

Prolonged use may cause headache and nausea(a type of Healing reaction). Use the product according to your health condition.

Acupressure marks are naturally left on the ears after use.

Do not use after drinking alcohol or when you are not feeling well.

When charging this product, be sure to use the dedicated charger that comes in the product box.

Sudden relaxation under excessive tension or stress may cause symptoms of autonomic dysfunction. Refrain from using this product under excessive tension or stress.

Do not throw or apply pressure to this product. Do not keep this product in a humid place, near a fire, or in a hot place.

This product is for home use. Do not use it for commercial purposes.

Do not wash the product with water. When cleaning the product, wipe it with a soft cotton cloth or with sterile cotton soaked in neutral detergent.

Do not use this product for treatment or diagnosis.

When this product is used by preschoolers or children, this product should be used under the guidance of parent or adult.

If you are under medical treatment or need medical advice, including those specified below, be sure to consult with doctor before use or refrain from using this product.

•Patients with restricted mobility or physical or mental disabilities
•Those with artificial heart pacemaker or various electronic aids
•Pregnant women
•Those receiving treatment of cervical disc disease

•Patients with serious diseases such as cancer, brain disease, heart disease
•If you are being treated for ear disease, or have any wound, inflammation or skin disease in the ears
•Those who are on medication such as taking high blood pressure medicine
•Those diagnosed with autonomic dysfunction

Healing Reaction

Healing reaction after using IAN is the evidence of improvement.

Healing reaction is a temporary worsening of symptoms that occur when the body is going through the process of healing itself through the process of finding the balance of the body. When symptoms improve, it disappears soon. The reaction is as follows. The timing and degree of the reaction depend on the individual.

01. Fever in ears.
02. Acupressure point is sore.
03. Ear is itchy.
04. Sleepiness.
05. Suddenly gets hot.
06. Headache, weariness, or body aches.
07. Feel languid or tired.
08. Sore throat, runny nose, or sputum.
09. Yawning, farts, or trim come out.
10. Temporarily nausea, vomiting, or diarrhea.
11. The amount of urine increases or the urine smells.

Troubleshooting Methods

Symptoms	Expected Causes	Troubleshooting Methods
No power	Insufficient charge or discharge	Charge the battery, or use the product after connecting the dedicated charger.
Weak acupressure	Insufficient charge	Charge the battery, or use the product after connecting the dedicated charger.
Not working	Bad connection of dedicated charger	Check to see if the dedicated charger is connected properly.

It can be used after 2 ~ 3 minutes when the dedicated charger is connected after full discharge.

※ For other inquiries, please contact our main contact number(82-1644-7671).

Product Specification

Model	IE-1000		
Rated Voltage	DC 5V		
Dimension	Size: 250(L) X 300(H) X 110(W) mm / Weight : 1,120g		
Optimal Working Environment	Temperature : 0°C ~ 35°C / Humidity: 80% or less		
Power Consumption	15W	Material	PC, ABS, Synthetic leather
Bluetooth	4.0	Battery Type	Lithium-ion battery

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